# Oclean

## Product User Manual

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## FCC WARNING

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

- 1.This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- 2.Children should be supervised to ensure that they do not play with the appliance.
- 3.It must only be supplied at SELV.
- 4.This appliance contains batteries that are non-replaceable.
- 5.The product can only be charged with the standard charging base. We will not responsible for any consequences resulting from the use of non-certified thirty charging base.
- 6.This appliance is not a fixed appliance, and the adhesive of the base is not used as a fixed appliance, but only as an auxiliary stabilizing effect.

EU REP

EU Authorized Representative  
Agency name: Mbrush UAB  
Representative: Algirdas Tevelis

European Address: Totoriug.3, Kaunas,  
Lithuania  
Post Code: LT44237

Email: [service@mobal.lt](mailto:service@mobal.lt)  
Phone: +37067052772

# Guarantee Certificate of Oclean Air2T Sonic Electric Toothbrush

## User Information

Name: \_\_\_\_\_ Phone: \_\_\_\_\_ E-mail: \_\_\_\_\_ ZIP code: \_\_\_\_\_

Address: \_\_\_\_\_

## Retailer Information

Retailer: \_\_\_\_\_ Phone: \_\_\_\_\_ Address: \_\_\_\_\_

Sales date: \_\_\_\_\_ ZIP code: \_\_\_\_\_ Invoice no.: \_\_\_\_\_

## Product Information

Product type: \_\_\_\_\_ Manufacturing site: \_\_\_\_\_

## Warranty Card of Oclean Air2T Sonic Electric Toothbrush

| Delivery date | Fault description/reason/disposal | Equipment employed | Acceptance date | Serviceman |
|---------------|-----------------------------------|--------------------|-----------------|------------|
|               |                                   |                    |                 |            |
|               |                                   |                    |                 |            |

Please request a maintenance report from the maintenance center and keep it together with the warranty certificate. (See the page on the left.)

Shenzhen Yunding Information Technology Co., Ltd., in strict accordance with the Responsibility Regulations on Maintenance, Replacement and Sales Return of Household Appliances, offers a guarantee for consumers who purchase the Oclean Smart Sonic Electric Toothbrush. In the case of any inconsistency between the service articles and the provisions of guarantee, the provisions of the guarantee and the contents of the guarantee certificate shall govern.

If the product you have purchased is found to have a fault or defect, please contact our service line at 4008-365-909.

Manufacturer: Shenzhen Yunding Information Technology Co., Ltd.

Manufacturer's address: 28G, Building 3, Dachong Business Center (phase III), No.18 Dachong 1st Road, Dachong Community, Yuehai Street, Nanshan District, Shenzhen, Guangdong, China.

Production address: No.89, Dingshan Middle Road, Dongfu Street, Haicang District, Xiamen, Fujian, China

Website: www.oclean.com

| Name of parts  | Toxic and harmful substance or element |    |    |      |      |       |              | Environment Protection |
|----------------|--|----|----|------|------|-------|--------------|------------------------|
|                | Pb                                     | Hg | Cd | Cr6+ | PBBs | PBDEs | 4 Phthalates |                        |
| PCBA           | ○                                      | ○  | ○  | ○    | ○    | ○     | ○            | 10 years               |
| Plastic parts  | ○                                      | ○  | ○  | ○    | ○    | ○     | ○            |                        |
| Metal hardware | ○                                      | ○  | ○  | ○    | ○    | ○     | ○            | 5 years                |

| Name of parts   | Toxic and harmful substance or element |    |    |      |      |       |              | Environment Protection |
|-----------------|--|----|----|------|------|-------|--------------|------------------------|
|                 | Pb                                     | Hg | Cd | Cr6+ | PBBs | PBDEs | 4 Phthalates |                        |
| Vibrating motor | ○                                      | ○  | ○  | ○    | ○    | ○     | ○            | 10 years               |
| Accessories     | ○                                      | ○  | ○  | ○    | ○    | ○     | ○            |                        |
| Battery         | ○                                      | ○  | ○  | --   | --   | --    | --           | 5 years                |

## Remarks

- This indicates that the content of toxic and hazardous substances in all homogeneous materials of this component is less than the limited quantity as specified in RoHS 2.0 Directive: 2011/65/EU & (EU) 2015/863.
- × This indicates that the content of toxic and hazardous substances in at least one of the homogeneous materials of this component is more than the limited quantity as specified in RoHS 2.0 Directive: 2011/65/EU & (EU) 2015/863.
- Refer to battery Directive: 2006 / 66 / EC and amendment 2013 / 56 / EU
- Except for very few materials, due to technical development reasons, most of the materials used are environmentally friendly. Please use this product with confidence during the environmental protection period.

FDA Registration Number: 3014246581  
CAN ICES-003(B)/NMB-003 (B)

Made in China